

## LATEX SURGICAL GLOVES POWDERED

12.0g

SLP 120

### Product Features

### Latex Surgical Gloves

Single Use	Pre-Powdered
Sterile	Natural Rubber Latex (NR)
Disposable	Hand Specific
	Curved Fingers
	Beaded Cuff
	Palm and Fingertip Textured
	Complies with USP Absorbable Dusting - Powder

Color Natural White, Off White

Size 7.5

also available in :  
5.5, 6.0, 6.5, 7.0, 8.0, 8.5, 9.0, 9.5, 10.0

### Physical Dimensions

Width (mm)	95 ± 4
Weight (g)	12.0g ± 1.0
Length	Min. 300mm / 12"
Mils	8
Palm Thickness (mm)	0.22 ± 0.02
Finger Thickness (mm)	0.23 ± 0.03
Cuff Thickness (mm)	0.14 ± 0.02

### ASTM Standard

Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	24 - 28	18 - 26
Elongation (%)	750 - 800	560 - 750
Modulus at 500% (MPa)	4.0 - 5.5	4.0 - 5.5

### Inspection Levels & AQL

### Inspection Level

### AQL

Dimensions	S2	4.0
Physical Properties	S2	4.0
Airblow Test	G1	0.65/1.0/1.5
Visual Inspection (Major)	G1	2.5
Visual Inspection (Minor)	G1	4.0

Protein Content 100 - 160 ug/dm<sup>2</sup>

Powder Content Below 12mg/glove

Major Certificates FDA

ISO 13485

MDA CLASS A

Country Specific Certificates

Test Compliance Reports Pre-shipment Reports

### Packaging

Box Pack (Sterile)

Bulk Packed (Non-Sterile)

### Storage Conditions

Store in a cool dry area with  
Temperature below 40°C

### Shelf Life

Product shall have a shelf life of  
3 years from the date of  
manufacture with the above  
storage conditions

### EN Standard

Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	12.0 - 14.0	9.0 - 11.0
Elongation at Break (%)	N/A	N/A

### Danars Healthcare Sdn. Bhd.

No. 27 Road 2/1,  
Tropicana Golf and Country Resort,  
47410,  
Petaling Jaya, Selangor,  
Malaysia.

W : [www.danarshealthcare.com](http://www.danarshealthcare.com)

T : +60 12 543 0740

E : [danarshc@gmail.com](mailto:danarshc@gmail.com)

SLP 120  
F8 1/1

25/04/23

