

Product Features

Cleanroom Nitrile Gloves

Single Use	Acrylonitrile Butadiene Copolymer Latex
Non- Sterile	High Titanium Dioxide (TiO ₂) dosage
Disposable	Powder Free
	Fingertip Textured
	Beaded Cuff
	Single Chlorinated
Color	White - Translucent
Size	XS , S , M , L , XL

Physical Dimensions

Width (mm)	XS	75 ± 5
	S	85 ± 5
	M	95 ± 5
	L	105 ± 5
	XL	115 ± 5
Weight (g)	XS	2.9g ± 0.2
	S	3.2g ± 0.2
	M	3.5g ± 0.2
	L	3.9g ± 0.2
	XL	4.2g ± 0.2
Length	All Sizes	Min. 240mm / 9"
Palm Thickness (mm)	0.07 ± 0.02	
Finger Thickness (mm)	0.10 ± 0.02	
Cuff Thickness (mm)	0.05 Min.	

ASTM Standard

Physical Properties	Before Ageing	After Ageing
Tensile Strength (Mpa)	Min. 14	Min. 14
Elongation (%)	500	400

Inspection Levels & AQL

Inspection Level

AQL

EN Standard		
Dimensions	S2	4.0
Physical Properties	S2	4.0
Airblow Test	G1	1.5
Visual Inspection (Major)	G1	2.5
Protein Content	NA	
Powder Content	Max 2mg/glove	
Major Certificates	Country Specific Certificates	
Test Compliance Reports	EN455-2	
	ASTM D6319	
	ASTM D6124	
	Trace Elements Test Report	
	Surface Resistivity Test Report	
	Static Decay Test Report	
	Pre-Shipment Reports	

Packaging	Box Packed
	Bulk Packed
Storage Conditions	Store in a cool dry area with Temperature below 40°C
Shelf Life	Product shall have a shelf life of 3 years from the date of manufacture with the above storage conditions

EN Standard

Physical Properties	Before Ageing	After Ageing
Force at Break (N)	Min. 6	Min. 6



Danars Healthcare Sdn. Bhd.

No. 27 Road 2/1,
Tropicana Golf and Country Resort,
47410,
Petaling Jaya, Selangor,
Malaysia.

W : www.danarshealthcare.com

T : +60 12 543 0740

E : danarshc@gmail.com